
# **Determining the Field of View, Size of Specimen & Magnification**

Field of View (FOV): the area you observe when you look through the microscope.

- Field diameter (FD): the diameter of the field of view
  - Often FOV and FD are used interchangeably.
- **NOTE –** As the magnification increases, the field diameter decreases proportionally.

As magnification **increases**, field of view **decreases** by the same factor.



As magnification **decreases**, field of view **increases** by the same factor.

		-	~ ~	
1mm	-	10	nn	um

Name of Lens	Objective Magnification (X)	Ocular Magnification (X)	*Total Magnification (X)	Field Diameter	
				mm	μm
Low Power					
Medium Power					
High Power					
	<b>*Total Mag</b> = Objective Mag x Ocular Mag				

# **Determining the Field of View (Field Diameter):**

- 1. Put the low-power objective lens in place. Place a transparent ruler on the stage. Position the millimetre marks of the ruler immediately below the objective lens.
- 2. Use the coarse-adjustment knob to focus the marks of the ruler.
- 3. Move the ruler so that one of the millimetre marks is just at the edge of the FOV. Note the diameter of the FOV, in millimetres, under the low-power objective lens. Then convert to micrometres (µm)
- 4. Rotate to the medium-power objective lens. Repeat steps 2 & 3 to measure the FD for this field. Make sure to convert to micrometres (μm)
- 5. Most high power objective lenses provide a field of view that is less than 1 mm in diameter, so it cannot be measured with a ruler. You must use the formula below.

#### $\rightarrow$ Calculating the Field of View for High-Power Lens

Calculate the ratio of the magnification of the high-power objective lens to the magnification of the low-power objective lens

Ratio = <u>Magnification of high-power lens (X)</u> Magnification of low-power lens (X)

Then use the ratio to determine the diameter of the field under high power magnification

**FOV (high power)** = <u>FOV low power (µm)</u> Ratio

### **Estimating Size of Specimen**

- 1. Choose the magnification that gives you the clearest image
- 2. Note the FOV/FD associated with that magnification
- 3. Estimate the number of times the specimen could fit across the field
- 4. Estimate the size of the specimen using the formula:

Estimated size of specimen = <u>Width of field of view (micrometres)</u> Number of specimens that fit across the field

\* Actual size = <u>field diameter (micrometres)</u> # of cells / specimens

### **Drawing Magnification**

The biological drawing magnification tells you how the size of the illustration compares with the actual size of the object.

Drawing Magnification = <u>Measured size of drawing (width or length?) (mm)</u> Estimated size of specimen (width or length?) (mm)


## **Practice Calculations**

- 1. Convert the following measurements:
  - a) 789.35 µm = \_\_\_\_\_ mm
  - b) 0.645 cm = \_\_\_\_\_ μm
  - c) 51.23 mm = \_\_\_\_\_ μm
  - d) 3.78 µm = \_\_\_\_\_ cm
- 2. Ben observed an *amoeba* under the 40X objective lens and noted that it took up 1/4 of the field of view. Calculate the actual size of the cell (in  $\mu$ m).

- 3. Ava examined a sample of pond water and found several protists in the water. She observed that 5 *Paramecium caudatum* cells could fit across the field diameter (FD) on high power. She followed the instructions and completed a formal drawing of one *P. caudatum* cell. The drawing measured 10.3 cm in length.
  - a) Calculate the estimated length of one *P. caudatum* cell (in  $\mu$ m)

b) Calculate the drawn magnification of the formal drawing.

- 4. A student observes an *Amoeba proteus* using high power magnification. The following diagram shows her observations:
  - a) Calculate the length of the Amoeba proteus in  $\mu m.$



b) The student then draws a diagram of the *amoeba* in her lab book. She uses a ruler to measure the length of her drawing, and determines that the drawing is 8.2 cm long. Calculate the drawn magnification of her *amoeba* drawing.